



Release 12 - Getting Started



Terms of use, trademarks, and copyright

All information in this document, including URL and other Internet Web site references, is subject to change without notice. Complying with all applicable copyright laws is the responsibility of the user. Without limiting the rights under copyright, no part of this document may be reproduced, stored in or introduced into a retrieval system, or transmitted in any form or by any means (electronic, mechanical, photocopying, recording, or otherwise), or for any purpose, except as provided in your license agreement, without the express written permission of Data Innovations LLC.

Data Innovations LLC may have patents, patent applications, trademarks, copyrights, or other intellectual property rights covering subject matter in this document. Except as expressly provided in any written license agreement from Data Innovations LLC, the furnishing of this document does not give you any license to these patents, trademarks, copyrights, or other intellectual property.

This documentation contains CONFIDENTIAL TRADE SECRET INFORMATION owned by Data Innovations LLC and/or its affiliated companies. Data Innovations and EP Evaluator are trademarks or registered trademarks of Data Innovations LLC.

Microsoft and Windows are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries.

All other company or product names are either trademarks or registered trademarks of their respective owners.

©1989-2017 Data Innovations LLC. All rights reserved.

Data Innovations LLC
120 Kimball Ave., Suite 100
South Burlington, VT USA 05403

Visit us at www.datainnovations.com

Welcome to EP Evaluator

Thank you for choosing this Data Innovations LLC product. This guide will lead you through the installation process for EP Evaluator.

Your EP Evaluator subscription includes:

1. Use of the software for 1-year.
2. Telephone support for the duration of the subscription.

North America: Monday through Friday, 9am-8pm EST/EDT:
(802)-658-2850

Latin America customer: +55-11-3801-3283

Europe: +32-2-332-24-13

Asia-Pacific: +852-2398-3182

3. Software updates (for the plan to which you subscribe) released during the subscription period.
4. Select Webcasts for training, please register for these trainings.

If you update EP Evaluator to the newest release, please notify Customer Support so that our records are accurate and you receive the correct activation information.

Getting Started

EP Evaluator Terminology

Major Release: A major release increments the release number. For example: from Release 10.0 to Release 11.0. A major release often addresses known issues and adds major enhancements or new functionality.

Minor Release: A minor release increments the build number. For example: from Release 11.1 to Release 11.2. A minor release often addresses known issues and adds minor enhancements.

Version: The version defines which of the two different types of EP Evaluator available for purchase is being referred to. The functionality of EP Evaluator is different between the **Professional** and **Standard** versions.

Activation Code: The activation code is provided by Data Innovations Customer Support by email when you purchase EE. The software is activated for the duration of your subscription.

Bring data forward: The process of restructuring EP Evaluator (EE) data gathered in an earlier major release to ensure it is compatible with the current release.

Software and Hardware Requirements

An IBM-compatible PC with a Pentium CPU or above.

For a 32-bit Microsoft Windows operating system:

- Windows 7

For a 64-bit Microsoft Windows operating system:

- Windows 7
- Windows 8.1
- 2008 Server
- 2012 Server
- Windows 10

NOTE: You will also require a PDF viewer to view EE reports.

Additional Specifications

- A minimum of 128 MB of RAM.
- Screen resolution: 1024 x 768 or greater.
- A minimum of 200 MB of hard drive space.

NOTE: The hard disk requirements of EP Evaluator as initially installed are less than 200 MB. If you generate lots of data (hundreds of projects), the total size required can exceed 2 GB.

NOTE: With EE Professional Version, the amount of disk space required may increase significantly if you use the Audit Trail feature.

Getting Started

Internet Access

Internet access is not required to operate EE. However, Internet access enables you to check for and download the latest EE build for your major release.

Locations EE Writes

EE users must have full read/write/delete access to any folder it writes to. Once installed and activated, EE only writes to three places on the computer:

- The directory tree where EE is installed
- The Application Data folder for the EE user (DGRhoad-sTempFiles)
- The HKEY_CURRENT_USER area of the registry

With the Professional Version of EE, an administrator can create data folders in alternate locations, including network locations. In this instance, EE can write to these additional locations.

NOTE: Every time EE runs, it performs a test to see if it can create and access subfolders and the Registry. If the test fails, EE displays an error message.

Note for Windows 7, 8.1 and 10 Users

If you decide to install EE in Program Files or Program Files(x86), Windows may write EE data to the virtual store area of the Windows Operating System. Although not generally a problem, this can lead to problems finding the data when using the Zipper program, after upgrading to a new major release of the software, or when browsing to the data using Windows Explorer.

Getting Started

Install EP Evaluator on a PC



NOTE: If installing on a network, see the following section on network installation before running the install wizard.

Access the EE setup file from your EP Evaluator CD-ROM or by downloading the software from the Data Innovations website.

If you already have an activation code use this link to download:

www.datainnovations.com/eedownload

Otherwise select the trial software download link on the following page and submit the request form. You will be sent a temporary activation code.

www.datainnovations.com/ep-evaluator

NOTE: The person installing EE must have administrator rights and must instal EE as administrator.

NOTE: If installing on a network, see the section on network installation before running the install wizard.

Install Process

INSTALLATION NOTE: For a minor release within major release, always install a later build on top of an earlier build. For example, install EE11.2 in the same folder as EE11.0 or EE 11.1. Never install a major release in the same directory as another major


release. For example, never install EE12 in the same folder as EE10, EE9, etc.

1. Right-click on the set up file with the name similar to Setup EE12-0(35).EXE.

NOTE: The filename of the EXE will contain different build numbers but will always use this basic format: SetupEEXX-d(dd).exe, where XX is the major release and d is a digit.

2. Select **Run as administrator** from the right-click menu.
3. The Setup Wizard launches. Click **Next** to start the installation process.
4. From the License Agreement window, use the scroll bar to review the terms of the End User License Agreement.
5. Select **I accept the terms of the license agreement** radio button and click **Next** to confirm your agreement.
6. From the Select Destination Location window, click **Next** to accept the default installation path of C:\EEXX (recommended).

Getting Started



7. Click **Finish**. EE is now installed on your machine and a desktop icon and link will have been created.

To start EE and run as Administrator

It is best practice to run EE as an administrator. You can configure compatibility properties on the program to always run as an administrator.

To configure compatibility properties:

1. Right-click the EP Evaluator icon on your desktop.
2. Select **Properties** from the right-click menu.
3. From the Compatibility menu, open the **Privilege Level** options.
4. Check **Run this program as an administrator**.

5. Click **Apply**.
6. Click **OK**.


You have configured the operating system to run EE as an administrator.

NOTE: If multiple users access EE from the same computer, update compatibility properties for all applicable user profiles.

Getting Started

Now Start EE

1. Double click on the desktop item to start EE.
2. If you have an activation code click **Activation Code**. If you want to learn how to obtain an activation code click **How to Order**.



3. The activation form will now be displayed. You must enter your organizational details on this page.

NOTE: These details cannot be changed once they have been entered and they will appear on project reports. Please ensure that accurate details are included here.

Enter Unlock Code X

| | |
|-------------|---------------------------------|
| Department | Hematology |
| Institution | Los Angeles University Hospital |
| Location: | Center 3H |

Activation Code: **SP02-2222F222-10APR2017-AAAA-BBBB**

Please fill in the Department, Institution, and Location information, which appears on all EP Evaluator reports. This is the only chance you will have to enter this information.

4. Click OK to start EE.

Getting Started

Network Installation

If installing EE on a network, you must create a **stable** shortcut to the software on the server, from each client machine.

NOTE: You must use the UNC (Universal Naming Convention) path for the shortcut rather than a mapped drive letter. We advise using the same shortcut on the server and clients. Right click on the desktop icon, select properties and type:

`\\servername\EE12\EE12.exe`

DO NOT USE, for example:

`Y:\EE12.exe`

For a network installation, use the same path for the server installation and for the shortcuts to run EE from the clients on the network.

Activation code for network installation:

When requesting an activation code for a network installation, the eight-digit System ID is different for each client. The first four digits identify the server and the last four digits uniquely identify the client (workstation).

For a network license, a client can access EE on a server as long as the first four digits of the client eight-digit System ID match the first four digits of the 8-digit System ID used to generate the activation code (server). If the first four digits of the server and client system IDs do not match, request an activation code for what's reported on

the client machine, and enter the activation code at that client. EE will work from the client, but not the server.

NOTE: For a network installation, it is necessary to activate EE as administrator from one client (workstation). It is also necessary to maintain a stable network path. If the location of the software on the network changes, the software will have to be re-activated for the new location and for each client.

Getting Started

Anti-Virus Software

If your company uses virus protection software, it is important to exclude the EE installation folder (EE11, EE12 etc.), and the **DGRhoadsTempFiles** directory from anti-virus scans. During a scan, anti-virus software may “lock” and/or modify the files being scanned. EE accesses the temporary files directory to store the temporary files necessary for the creation of reports and print previews. The location of the DGRhoadsTempFiles directory is dependent upon the operating system:

Windows 7:

C:\Users\

Windows 8.1

C:\Users\

Windows 2008 Server

C:\Users\

Windows 2012 Server

C:\Users\

Windows 10

C:\Users\

NOTE: The location of the folder is dependent on your system. The files may be hidden on your system. To show hidden files: Open Folder Options by clicking the Start button. Click Control Panel, Appearance and Personalization, Folder Option, View tab. Under Advanced settings, click Show hidden files and folders, and then click OK.

Activate a trial version

EP Evaluator comes with a free 30-day trial period. Once this trial period passes, you will be prompted to activate the software.

Data Innovations will provide you with an Activation Code when you purchase the software.

Follow the steps below to activate the software:

1. Launch EP Evaluator.
2. Click the Activate button on the Notice screen.
3. Update the **Department**, **Institution**, and **Location** details.
4. Enter the Activation code in the **Activation Code** text box.
5. Click **OK**.

The software will confirm that you have activated EE.

Upgrading EP Evaluator

Minor Release

For a minor release within a major release, always install a later build on top of an earlier build. For example:

- Install EE11.2 in the same folder as EE11.0 or EE 11.1.

For minor releases, the activation code does not change and all data is automatically restructured to ensure consistency with the more recent build.

Major Release

If moving to a new major release, never install the new major release in the same folder as an older major release. For example:

- Never install EE11 in the same folder as EE10, EE 10.1, etc.

When moving to another major release of EE, you must contact Data Innovations support to purchase a new Activation Code.

Data will need to be manually brought forward from the older major release into the newer major release.

NOTE: EE Professional version users should contact support before installing a newer release.

Bring data forward to a new major release

Bringing data forward manually ensures that data from the older release is structured appropriately for the newer release. There are two options for bringing data forward:

1. The Upgrade Wizard option brings data forward from all active studies in the earlier release.
2. The Backup and Restore option allows you to select which backed-up (archived) projects from the earlier release should be restored in the new release.


To bring data forward from one computer to another, you must use the Backup and Restore option to backup data on prior releases before restoring them on the current EE install.

NOTE: Backing up completed projects on a regular basis streamlines the process of bringing data forward using the Backup and Restore option, because most of the data you wish to bring forward is already archived.

Upgrade Wizard

1. Ensure the data from the prior EE release is on the PC where the new release is installed.
2. Launch the new release of EE.
3. From the **Statistical Module** window, open the **Utilities** menu.
4. Click on **Update Wizard**.
5. From the Update Wizard window, select the EE program you are moving data from.

Getting Started



NOTE: To bring data forward from any release prior to EE4, using the update wizard, contact Data Innovations support.

6. Browse to the STUDIES folder that contains the data.

NOTE: The final folder in the directory path must be called STUDIES (not case-sensitive).


7. Click Finish.

NOTE: The Update Wizard will move all your data forward into the new EE release. If you are bringing forward a large amount of data, this process can take a significant amount of time.

A dialog confirmation screen indicates the Update Wizard has successfully completed bringing data forward. The Activity Log displays additional information or warnings.

Backup and Restore

1. Launch the older release of EE.
2. Open the **Utilities** menu.
3. Open **File Manager**.



4. Select the project(s) you wish to back up. To back up a group of projects, select the group prefix in the Project list. All projects in the highlighted group will be backed up. To select several projects at once hold down the shift key before clicking with the mouse pointer. To choose individual projects hold down the control key before selecting the projects.

5. Click **Backup**.
6. Repeat steps 2 through 5 for all projects you wish to back up.
7. Copy all backed-up projects from the \DATA\Backups folder for the older release to the \DATA\Backups folder for the new release.
8. Close the older release of EE.
9. Launch the new major release of EE.
10. Open the **Utilities** menu.
11. Open **File Manager**.
12. From the **Backup Files** box, select a backup file to restore.
13. Click **Restore**.
14. Repeat steps 10 through 13 for each project you wish to restore to the new major release

Getting Started

User Help

Electronic help can be accessed from within EP Evaluator by pressing F1 or from the EE Help menu.

Additional documentation can be downloaded from our Website:

<http://www.datainnovations.com/ep-evaluator>

DATA INNOVATIONS NORTH AMERICA

120 Kimball Ave., Suite 100

South Burlington, VT USA 05403 Fax: (802) 658-2782

Support: Phone: (802) 658-1955

Email: northamerica-support@datainnovations.com

Sales: Phone: (802) 264-3470

Email: northamerica-sales@datainnovations.com

DATA INNOVATIONS ASIA

Fax: +852 2398 8667

Support: Phone: +852 2398 3182

Email: asia-support@datainnovations.com

Sales: Phone: +852 2398 3182

Email: asia-sales@datainnovations.com

DATA INNOVATIONS LATIN AMERICA

Fax: +55 (11) 3871-9592

Support: Phone: +55 (11) 3801-3283

Email: latinamerica-support@datainnovations.com

Sales: Phone: +55 (11) 3801-3283

Email: latinamerica-sales@datainnovations.com

DATA INNOVATIONS EUROPE

Fax: +32 2 376 43 84

Support: Phone: +32 2 332 24 13

Email: europe-support@datainnovations.com

Sales: Phone: +32 2 770 62 22

Email: europe-sales@datainnovations.com

